


Adapter Kit for ADEC® Bottle System



Easy mounting to an existing ADEC® bottle system

- Water bottle receptacle for ADEC® bottle system
 - easy conversion of the ADEC® bottle system to the Medtronic bottle system
 - water bottle is connected to compressed air (regulated by an automatic pressure regulator)
 - air pressure stops automatically when water bottle is exchanged



Medtronic water bottle receptacle for ADEC® bottle system

One contact only!

We have our own development department, manufacture our own tools and moulds and produce over 90 % of our products ourselves. Therefore, we are in a position to respond to any new product requirements quickly and can ensure highest quality standards. Medtronic is certified according to DIN EN ISO 13485.