


MC 3FP/3FP-S Syringe Angled or Straight Form



! Advantages

- Very light! Weighs only 35g (without hose)
- Syringe head can be removed for easy cleaning/disinfection
- with soft grip handle to avoid slipping during use
- Quick-disconnect feature
- To be used with disposable syringe tips, Stropko Irrigator, autoclavable metal syringe tips
- CE approved

Soft Grips available in the following colours:



For further details please see reverse!


3-Way Syringe

Details

- Weighs **only 35 g!**
About 2/3 less than regular metal syringes!
- Quick disconnect feature for mounting disposable syringe tips, Stropko Irrigator, autoclavable metal syringe tips
- Clear separation of air and water, according to European standard (air on the left; water on the right) or US standard (air on the right; water on the left)
- Push buttons for proper dosing of air and water
- Almost no release of heat due to use of HTC polyacetal resin
- Soft grip to avoid slipping and ensure safe use
- Service-friendly control valves
- No metal abrasion on instrument holder
- Company logo may be imprinted (upon request)
- **Optional:** Syringe head removable for easy cleaning; Quick shut-off for water and air
- CE approved



! Sleeves available for different types of instrument holders (please inform us which syringe has been used previously)!



One contact only!

We have our own development department, manufacture our own tools and moulds and produce over 90 % of our products ourselves. Therefore, we are in a position to respond to any new product requirements quickly and can ensure highest quality standards. Medtronic is certified according to DIN EN ISO 13485.

medizinisch-elektronische Geräte-Gesellschaft m.b.H.

Riedelbach | Technologiepark
In den Ensterwiesen 7
61276 Weilrod | Germany

Tel +49 (0) 6083 95989-0
Fax +49 (0) 6083 95989-11
info@medtronic-dental.de
www.medtronic-dental.de