


## MedtroSmart Dental Unit for Air Driven Instruments

- Small, robust and low-maintenance dental unit
- Includes:
  - 1 x MC 3 FP-S syringe
  - 2 x 4 hole turbine connection (for turbines, air motors or air scaler)
  - instrument hoses including oil collector and unit side connector
  - separate water bottle system „Clean Water“
  - pneumatic foot switch with unit side connector
- Anodised aluminium housing
- With wall fastening or mounting bracket
- For use only air pressure necessary
- Certified medical device
- Company logo may be laser engraved (upon request)

*(Instruments not included in delivery)*




### One contact only!

We have our own development department, manufacture our own tools and moulds and produce over 90 % of our products ourselves. Therefore, we are in a position to respond to any new product requirements quickly and can ensure highest quality standards. Medtronic is certified according to DIN EN ISO 13485.

## MedtroSmart Plus Dental Unit for Air Driven Instruments

- Small, robust and low-maintenance dental unit
- Includes:
  - 1 x MC 3 FP-S syringe
  - 2 x 4 hole turbine connection (for turbines, air motors or air scaler)
  - instrument hoses including oil collector and unit side connector
  - separate water bottle system „Clean Water“
  - pneumatic foot switch with unit side connector
- Corian housing
- With wall fastening or mounting bracket
- For use only air pressure necessary
- Certified medical device

*(Instruments not included in delivery)*



### One contact only!

We have our own development department, manufacture our own tools and moulds and produce over 90 % of our products ourselves. Therefore, we are in a position to respond to any new product requirements quickly and can ensure highest quality standards. Medtronic is certified according to DIN EN ISO 13485.